

Discovery of antigen-specific B-cell and T-cell clonotypes using a multiplexed dCODE Dextramer®-based workflow

Liselotte Brix¹, Kivin Jacobsen¹ and Dilek Inekci¹
¹Immudex, Copenhagen, Denmark

Introduction

Understanding the antigen-specific B and T cell responses is key for development of vaccines and targeted therapies, encompassing various stages from target discovery to monitoring the treatment efficacy to patient stratification. The DNA barcoded dCODE Dextramer® and dCODE Klickmer® reagents allow detection of antigen-specific cells through their corresponding B and T cell receptors (BCR/TCR), providing access to full-length paired receptor sequences, gene expression profile, and cell surface proteins from the same cell when used in combination with single cell technologies. Obtaining full-length, paired BCR sequences provides a resource for efficient antibody discovery for therapeutic and prophylactic use. Likewise, accessing full-length, paired TCR sequences enables investigation of receptor potential in adoptive cell therapies like CAR T-cell therapy.

Here, we established a workflow for simultaneous discovery of antigen-specific B-cell and T-cell clonotypes using a multiplexed dCODE Dextramer®-based approach combined with the 10x Chromium Single-Cell Analysis System.

Workflow for simultaneous investigation of antigen-specific B and T cells in PBMC samples






Figure 1. Experimental overview. A virus-specific model was used to establish the workflow containing a panel of dCODE MHC Dextramer® and dCODE Klickmer® reagents covering 21 specificities for detection of antigen-specific B and T cells, respectively, recognizing **CMV, EBV, flu, SARS-CoV-2** antigens and controls. The workflow was demonstrated by assessing the virus specific cells and their BCR/TCR clonotypes after a SARS-CoV-2 infection. PBMC samples were collected from a donor immediately after (0-day) and six weeks after (6-week) infection. Each of the PBMC samples was **(1)** stained with the panel of dCODE® reagents in the same tube and subjected to **(2)** 2-way sorting of antigen-specific cells by flow cytometry to sort Ag-specific B and T cells, respectively, **(3)** pooling of sorted Ag-specific B and T cells, **(4)** 10x Chromium Single-Cell Analysis System, **(5)** Library preparation and sequencing and **(6)** BCR/TCR identification.

Divergences in the BCR repertoire shortly following infection versus six weeks post-infection



Divergences in the TCR repertoire shortly following infection versus six weeks post-infection



Conclusion

- (1) We have demonstrated a multiplexed workflow for simultaneous detection of Ag-specific B and T cells in blood samples also providing information on their corresponding BCR and TCR clonotypes, gene expression and surface marker profiles
- (2) We included a total of 21 antigen specificities using a virus-based model system
- (3) The workflow can be tailored to cancer-specific responses by using dCODE® reagents displaying cancer antigens
- (4) The workflow allow characterization of disease-specific immune responses with potential to facilitate diagnostic approaches, prediction of disease progression and development of novel cancer vaccines and therapeutics