

Dissecting MAIT cell phenotype and TCR sequences using MR1 reagents and single-cell RNA sequencing

Bjarke Endel Hansen¹, Thomas Holberg Blicher¹, Kivin Jacobsen¹, Liselotte Brix¹


¹Immudex Aps, Virum, Denmark

Introduction

Mucosal-associated invariant T (MAIT) cells constitute a large and specialized subset of T cells with both innate and adaptive features. Unlike conventional T cells, MAIT cell TCRs exhibit reduced sequence variability and target metabolites of both bacterial and endogenous origin in the context of the MHC-I-like molecule known as MR1. To help study MAIT cells we have developed MAIT cells specific reagents, MR1 dCODE Dextramer®, and demonstrate here:

- a workflow, which enables the study of MAIT cell phenotypic markers combined with TCR sequencing using MR1 dCODE Dextramer® reagents on the BD Rhapsody™ Single-Cell Analysis System.
- that recombinantly expressed TCRs from identified MAIT TCR sequences specifically bind to their antigen, MR1/5-OP-RU, in a bead-based assay.

Workflow of the full immunological profiling of MAIT cells



The full immunological profiling of immune cells in a single workflow by combining dCODE Dextramer® reagents and the BD Rhapsody™ Single-Cell Analysis System.

Conclusion

- MR1 dCODE Dextramer® reagents achieve sensitive and specific detection of MAIT cells and, when used with the BD Rhapsody™ Single-Cell Analysis System, facilitate single-cell multi-dimensional analyses that reveal correlations between phenotype, TCR sequence, and gene expression.
- The MR1/5-OP-RU dCODE Dextramer® reagent reliably and specifically identified CD3+CD8+CD161^{high} that phenotypically were memory-like (CCR7-CD45RA^{low}CD62L^{low}).
- Differential gene expression of MR1/5-OP-RU-specific MAIT cells compared to other CD3+ T cells reflected their unique response to immunological environments and signals.
- The V(D)J and CDR3 analysis confirmed the reported semi-invariant TCR α chain of MAIT cells. TRAV1-2, TRAJ33, TRBV6, and TRBV20 were the most used gene segments. Despite both α and β chains of the TCR showed restricted gene segment usage, the TCR repertoire of the MAIT cells was diverse, with no major clonotypic dominance.
- Recombinant expression of two MAIT TCRs confirmed their specificity to MR1/5-OP-RU, but also indicated an antigen independent binding to MR1.

Accurate detection of MAIT cells within the sequencing data




Figure 1. UMAP clustering based on gene expression data. Overlaying clusters with surface markers (BD Rhapsody™ Immune response panel antibodies and MR1 dCODE Dextramer®) identified the major cell lineages (CD4+ and CD8+ T cells, B cells, Monocytes, NK cells) and CD161^{high} MAIT cells among the CD3+ T cells.




Figure 2. MAIT cells exhibit a memory-like phenotype (CCR7-CD45RA^{low}CD62L^{low}) and a non-active (CD137-CD25-HLA-DR-)(not shown) and response-ready state (CD27^{high}). **a)** MAIT cells were identified by cell surface markers (CD3+CD161^{high}MR1/Dex+). **b)** Cell surface marker analysis of MAIT cells.




Figure 3. Differential gene expression: shows the top up- and downregulated genes identified in the MAIT cell population compared to the remaining CD3+ T cells. Several of the differentially expressed genes are immune-related, indicative of the unique response of MAIT cells to different immunological environments and signals.

V(D)J analysis of MAIT TCR shows high variation




Figure 4. TRAV and TRBV gene segment representation in all TCR clonotypes (**A**) and the top 10 most frequent clonotypes (**B**) of the analyzed MAIT cell population. The majority of the identified CD161+ MAIT cells (~70%) express typical TCRs consisting of TRAV1-2 and TRAJ33(not shown) as well as TRBV6 and TRBV20 gene segments. Despite the restricted gene segment usage, overall paired TCR sequence variation was high, with 138 clonotypes in 158 cells. (**C**) gene segment usage and CDR3 of 10 most frequent clones.

Validation of MAIT TCR specificity in bead-based assay




Figure 5. Recombinant TCR proteins bind specifically to MR1/5-OP-RU.

Two identified MAIT TCRs (TCR032 and TCR033) were recombinantly expressed in E.coli, refolded and evaluated for binding to their target, MR1/5-OP-RU. Both MAIT derived TCRs, but not the control TCR (TCR002), bind efficiently to MR1/5-OP-RU conjugated beads, but not to control HLA-A*0201 conjugated beads. Some binding to the MR1 loaded with control metabolite 6-FP was also observed, and is probably due to general, but weak interaction with the MR1 protein independently of the bound ligand.